


FOR USE ON THE JAY GS® CUSHION

1. Remove the cover and fluid pad from the cushion to determine correct placement of the Lateral Thigh Support.
2. Detach and discard extra Velcro® strips from the Thigh Support. They are not required for use with the Jay GS Cushion.
3. Place Thigh Support on cushion base at desired location. Attach the large Velcro Suspenders (4 x 6") so that the suspender contacts both the 1 x 4" piece of loop Velcro on the Thigh Support surface and the loop Velcro on the cushion base. Press down firmly on the suspender.



FOR USE ON THE JAY® J2 OR JAY® CUSHION

1. Remove the cover and fluid pad from the cushion to determine correct placement of the Thigh Support.
2. Place Thigh Support at desired position on Cushion base. Using rubbing alcohol thoroughly clean cushion base where paper-backed Velcro pieces will be adhered. Remove the backing to the 1 x 4" loop Velcro piece and press Thigh Support back into position.
3. On the cushion base place the extra 1 x 5-1/2" loop Velcro piece parallel to and next to the Thigh Support edge.
4. Attach the large Velcro Suspenders (4 x 6") so that the suspender contacts both the 1 x 4" piece of loop Velcro on the Thigh Support surface and the 1 x 5-1/2" loop Velcro on the cushion base. Press down firmly on the suspender.

