


INSTALLATION AND ADJUSTMENT OF THE LATERAL THIGH SUPPORT

The left and right thigh support pads and brackets can be mounted to the asymmetrical leg plates, using 2x T-nuts, 2x washer, and 2x M6x1x12 bolts per side.

1. Prepare leg plates for installation

1 2

- You can use an utility knife, or peel away the excess velcro (A) in the case of a slot being covered up.
- Slots exposed, and velcro removed (B)
- It may be necessary to loosen the screws (C) under the leg plates so that they can slide forward to expose the outer slots

**2. Attach the lateral thigh support**

3 4 5

- The lateral thigh support and pad (D) can be installed in 3 different locations on the leg plate: Narrow, Medium, and Wide.
- Using a 4mm hex key, attach the thigh support bracket using 2x button head screws, 2x washers, and 2x T-nuts. Do not tighten the hardware until correct adjustment has been made.
- The bracket and pad can be rotated slightly to support most leg angles.
- Using the 4mm hex key, tighten the screws.
- If necessary, repeat the steps for the other side.

